



How to Package In-House Aligners

An AAO TechSelect Guide



1 Create label to identify each aligner

- Labels branded with your practice name are recommended to indicate patient name and aligner number in the series.



2 Package aligner in sealed bag.

- Depending on the bag size, package each set of aligners or each aligner separately in a sealed bag and label. Popular methods include cellophane bags that you heat seal or small Ziploc type bags.



3 Place packaged aligners in labeled box for chairside delivery

- Place aligner series in box labeled with patient information. This aids in easier transfer from lab to patient and from your office home with patient.
- Note: If you are using clear bags and UV light to sterilize, you can perform the sterilization step after the aligners have been packaged and labeled.

